

# Fast. Easy. Affordable.

Go Digital – with the Upgrade Solution that's Right for You.



SMART IMAGING SOLUTIONS

# Digital Retrofit Solutions for X-ray Rooms and Mobile Systems.

If your facility is still imaging with film or CR, an upgrade to full wireless digital can provide tremendous benefits: superb image quality, greater productivity and super-charged workflow. That's why Carestream offers a wide array of retrofit solutions – all designed to help you move your X-ray rooms and mobile units up to digital – quickly, easily and affordably. Our retrofit solutions allow you to leverage the investment in your existing equipment while gaining all the benefits of a wireless DR workflow. And you have a variety of detector options to choose from. Look to our Focus line of detectors for an affordable, ease-of-entry upgrade to full wireless DR. Or choose the CARESTREAM DRX Plus or Lux 35 Detectors that feature the X-Factor so your detectors can be shared with other DRX equipment in your facility.

## For Room Upgrades:

The CARESTREAM DRX-1 System offers an easy way to upgrade X-ray rooms to premium digital performance with a high-performing, wireless detector.

## For Room or Mobile Upgrades:

Here's our flagship upgrade. The CARESTREAM DRX-Transportable System/Lite delivers wireless tablet operation and superb digital image quality. And it's shareable across both your rooms and mobile systems.

## For Mobile Upgrades:

Upgrade mobile systems from nearly any manufacturer with the CARESTREAM DRX-Transportable System/Universal Mobile.

# Room Upgrades

## The CARESTREAM DRX-1 System.

Upgrade your X-ray room to DR in one day.

Converting to industry-leading DR quality, workflow and productivity is fast and easy with the DRX-1 System. The system is designed to help optimize the imaging equipment you're currently using.



Each detector simply slides into your existing X-ray equipment. The detectors meet ISO 4090 standards – so there's no need to modify your existing table or wallstand Bucky, replace your grids or radically alter your procedures. Plus, the DRX-1 System can be installed in your exam room in just one day.

- The detector wirelessly transmits digital images to your capture console.
- Images are ready for immediate viewing and manipulation.
- From the console, images can be quickly forwarded to a printer, PACS archive and other network destinations.
- The detector's lightweight, wireless operation provides convenient handling and excellent positioning flexibility.
- The advanced enclosure design protects against dust and water ingress with an IP57 rating for our DRX Plus and Lux 35 Detectors, and an IP44 rating for our Focus Detectors.
- State-of-the-art CARESTREAM ImageView Software, powered by Eclipse, delivers superb image quality and an enhanced workflow.



# Room/Mobile Upgrade

## The CARESTREAM DRX-Transportable System/Lite.

Our flagship DR room or mobile system retrofit solution.

A versatile, all-in-one solution to upgrade your rooms and mobile systems to full wireless DR.



You can use it in your room, then move it to your mobile system, then move it back to your room again. No need to invest in multiple systems. This makes upgrading to DR amazingly easy, and wireless tablet operation gives you unprecedented freedom. There's no bulky electronics box required, and no need to work from a fixed console. The lightweight, wireless tablet gives you complete control of the system and immediately displays your images after they're acquired.

- The Virtual Access Point in this detector makes it a flexible solution that enhances workflow.
- Quickly and cost-effectively convert existing analog rooms or mobile system to DR.
- Self-powered and rechargeable, its beam-sensing technology further simplifies your upgrade.
- This upgrade lets you move easily around the Radiology Department, from room to room or from room to mobile.



## Meet Our Family of DR Detectors.

Our portfolio of detectors includes wireless, shareable, glass-free and fixed models to meet your demands based on exam type, detector size, dosage level and budget.

Our DRX Plus and the Lux 35 Detectors feature the X-Factor, so they can be shared across all of your DRX products. So think sharing, think X-Factor to maximize your investment and improve productivity.

Or choose our Focus Detectors for superb image quality at a highly economical price.



Discover the right detector for your needs at [carestream.com/detectors](http://carestream.com/detectors).

# Mobile Upgrade

## The CARESTREAM DRX-Transportable System/Universal Mobile.

Designed to retrofit mobile systems from nearly any manufacturer.

The DRX-Transportable/Universal Mobile is a tremendous performance-booster.

- It offers a flexible design for mobile imaging systems from multiple suppliers.
- It's ideal for trauma sites, remote clinics, operating rooms and the ICU.
- Choose from multiple tablet PC mounting options to fit your mobile X-ray system.
- Electronics fit in the existing cassette bin.
- A 13-inch LCD touchscreen display provides 3000 x 2000 resolution.
- The system is self-powered and rechargeable.



NO MATTER WHICH ROAD YOU'RE ON,  
WE'LL GET YOU WHERE  
YOU NEED TO GO.

SMART IMAGING SOLUTIONS

### The Fastest Route to a Successful Imaging Future.

As you travel to the future of imaging, what's your next step? Perhaps you need to accelerate your X-ray workflow. Or add a new imaging solution to your facility to provide better patient care.

Wherever you're headed, we have the smart solutions you need: Full-digital imaging rooms and leading-edge mobile units. Wireless, shareable detectors and CR-to-DR retrofit kits. Affordable CR systems. Advanced digital imagers and our Managed Print Solutions service. And as always, the gold-standard in X-ray film technology.

### Start Mapping Your Route Today.

From small clinics, urgent care facilities and specialty practices to large hospitals and imaging centers, every facility will find a Carestream solution to meet the most specific needs. Visit [carestream.com](http://carestream.com) and start your journey to greater productivity and a higher standard of patient care.

Carestream also offers Non-Destructive Testing solutions to help ensure quality and safety in the aerospace, petrochemical and other industries, as well as Contract Manufacturing to bring innovative coated products to market.

### A Community of Service and Support.

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver.



### Eclipse: Harness the Power of Decades of Innovation.

Carestream's Eclipse is the engine behind our innovative imaging software. It uses AI technology and proprietary algorithms to significantly increase the value of the entire imaging chain, from capture to diagnosis.

[carestream.com](http://carestream.com)



"Rx only"

© Carestream Health, Inc., 2022. CARESTREAM is a trademark of Carestream Health. CAT 2000 273 02/22

Carestream